

- **The Shoulder Injection Trainer** For practicing soft tissue joint injection used for the treatment of injuries and arthritis. There are five different injection sites incorporated into this model
- Represents normal anatomy precisely (both superficially and subdermally) for palpation purposes
- Feedback Console illuminates:
 - green LED when needle is correctly placed in each of the injection sites
- Supplied with Joint and Soft Tissue Injection book by Trevor Silver
- Skin surface is washable using soap and water
- Injection on 5 specific sites:
 - subacromial space
 - acromioclavicular joint
 - bicipital groove
 - glenoid fossa from the anterior aspect
 - glenoid fossa from the posterior aspect



